

Cleaning and Infection Control of Devices in CF Patients

Eynav Manor, B.PT, M.A.

Head of Pediatric Respiratory Physiotherapy

Rambam CF center

Ruth Rappaport Children's Hospital

Rambam Health Care Campus



Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

General

- 📌 Cystic Fibrosis is a life shortening disease.
- 📌 There is decreased mucociliary clearance and impaired host defenses.
- 📌 Bacterial infection plays a major role in the process leading to respiratory failure in CF.

* Cohen-Cymberknoh et al, Am J Respir Crit Care Med 2011, Jun 183(11):1463-1471

* Ehre, Int J Biochem Cell Biol 2014, Jul 52:136-145

* Peckham et al, J Cyst Fibros 2015, Jun 20 pii:S1569-1993(15)00153-8



Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

General

✚ Respiratory physiotherapy:

- ✚ Improves and assists in airway clearance.
- ✚ Keeps airways hygiene.
- ✚ Improves medication delivery.

✚ Many of our CF patients use at least one respiratory device in the physiotherapy session.

* Rand, Paediatr Respir Rev. 2013 Dec;14(4):263-9



Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

General

- ✚ Most of the literature refers to the contamination and cleaning of nebulizers.
- ✚ It is unclear if a device can contaminate the patient.
- ✚ Staphylococci and enterococci survival on plastic more than 90 days.

* Peckham et al, J Cyst Fibros 2015 Jun 20;pii:S1569-1993(15)00153-8

* Saiman et al, Infect Control Hosp Epidemiol. 2014 Aug;35 Suppl 1:S1-67

* Saiman, Clin Microbiol Rev 2004 Jan;17(1):57-71

* Neely, J Clin Microbiol 2000 Feb;38(2):724-6

* Reyhner, J Cyst Fibros 2005 Sep;4(3):183-7



Organism transmission routes

Organism	Transmission routes
<i>Pseudomonas aeruginosa</i>	Person-to-person, Environmental reservoir, Environmental surface
<i>Burkholderia cepacia</i> complex	Person-to-person > Environmental surface
<i>Stenotrophomonas maitophilia</i>	Environmental surface, Environmental reservoir > person-to-person
<i>Achromobacter xylosoxidans</i>	Environmental surface > person-to-person
Non-tuberculous <i>Mycobacterium</i>	Environmental reservoir > > person-to-person
<i>Aspergillus</i> species	Environmental reservoir > > person-to-person + airborne

CFF guideline for nebulizers

- ✚ Cleaning nebulizer parts with dish detergent and water.
- ✚ Disinfect with one of the following options, if permitted by the manufacturer.
- ✚ Air-dry completely.

* O'Malley, Respir Care 2009 May;54(5):641-657

* Saiman et al, Infect Control Hosp Epidemiol. 2014 Aug;35 Suppl 1:S1-67



Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

Aim

- ✚ To determine if the respiratory devices can be contaminated following use.
- ✚ To determine if the cleaning guidelines for nebulizers are effective for respiratory devices.



Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

Method

- 📌 The patient brought his/her respiratory device to the clinic, following use without cleaning.
- 📌 Just the devices that can be taken apart.
- 📌 Since March 2015 and ongoing.



Method

- ✚ Swab (10 ml SRK swab, Copan Italy) before cleaning.
- ✚ Cleaning: The parts of the device were soaked in hot water (not boiling) + standard dish detergent, for 10 min
- ✚ Left to dry.
- ✚ Swab again from the same places.
- ✚ Sent to the laboratory.

* Saiman, Clin Microbiol Rev 2004 Jan;17(1):57-71



Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

Laboratory

- ✚ Swabs were transferred to the microbiology laboratory within 2 hours.
- ✚ 100 μ l aliquots were spread on Brain Heart Infusion agar plates and were incubated for 48h at 37°C with 5% CO₂.
- ✚ Total count CFU (colony forming units) was determined from the plates and predominant colonies were identified using MALDI-TOF MS technology (Vitek-MS, Biomerieux, France).



Ruth


Ruth Rappaport Children's Hospital

Rambam Health Care Campus


Method

Questionnaire:

 How often do you use your device?

 Once a day / several times a week / once a week / once in two weeks / once a month / not at all

 How often do you clean your device?

 Once a day / several times a week / once a week / once in two weeks / once a month / not at all

 How do you clean the device?

 Which cleaning agent do you use?



Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

Results

patient	Samples	Age Av.	PI	FEV1 Av.
18	20	16.64 (0.6-50)	13	63.47 (32-97%)

Devices	
IPV	11
Ressistex	3
Flutter	2
Aerobika	2
Pipep	1
Acapella choice	1



Results

Devices use	
Daily	11
Several times a week	6
Once in 2 weeks	2
Once a year	1

Devices cleaning	
Daily	8
Several times a week	2
Once a week	5
Once a month	2
Never	3




Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

Results



Ruth
Ruth Rappaport Children's Hospital
Rambam Health Care Campus

CFU After CFU Before

Results

- 📌 In 10 samples - cleaning completely eradicate all CFU.
- 📌 In 2 samples - no CFU before and after cleaning.
- 📌 In 6 samples - cleaning was partially effective.
- 📌 In 2 samples - cleaning was ineffective.

CFU before Av.	CFU after Av.	Paired T-test
63,586.5	50,965	0.036186



Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

Results

Bacterial Colonies

Before cleaning	After cleaning	Sputum culture
Staphylococcus pasteuri	ללא שינוי בכמות/ממצאים	Candida
Pseudomonas stuzeri, Staph warneri, Corynebacterium sp.	ללא שינוי באורגניזם ירידה של פי 100 בכמות המושבות	Candida, Achromobacter xyloso., staph aur.
Corynebacterium sp	ללא שינוי בכמות/ממצאים	Candida, Pseudomonas aer.
Moraxella sp	ללא צמיחה כלל	Candida, Aspergillus niger.
סטרפטוקוקים	ללא צמיחה כלל	Candida, Serratia marc.

Limitations

- ✚ Small sample.
- ✚ Cleaning and usage extent were reported by patients.
- ✚ The interval from the last cleaning was not determined.
- ✚ Environmental swabs do not enable diagnosis of specific bacteria.



Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

Conclusions

- ✚ Respiratory devices are contaminated following use.
- ✚ Appropriate cleaning reduces contamination.
- ✚ It is known that patients can be infected by contact with environmental surfaces (e.g pseudomonas, Burkholderia).
- ✚ It is plausible that patients can be cross infected from the devices.



Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

Conclusion

- ✚ The respiratory devices have to be cleaned after each use.
- ✚ There is a need for infection control guidelines of the growing list of respiratory devices.
- ✚ Recommendation - periodically refresh the device cleaning process.



Ruth

Ruth Rappaport Children's Hospital

Rambam Health Care Campus

Thanks

🏆 Prof. Lea Bentur

🏆 Dr. Michal Gur

🏆 Dr. Yuval Geffen

🏆 Team of Pediatric Pulmonary Division

🏆 All the patients

